



Professional disinfecting fluid · 5L

Milton Disinfecting Fluid

Sterilising solution for equipment & surfaces

15 min**No rinse needed****24 hours**

cold water sterilising on sterilised items solution stays effective

Milton Disinfecting Fluid is a sodium hypochlorite-based disinfectant for hospitals and professional kitchens — work surfaces, floors, walls and utensils. It kills bacteria, fungi and viruses. No heating and no rinsing needed once items come out of solution.

Before you start

- Do not use on metal — it can cause corrosion.
- Never mix with other products — this can release toxic chlorine gas.
- Avoid contact with eyes and skin; splashes may discolour fabric.
- Use a non-metallic container for the solution.
- Full safety information is in the product's Safety Data Sheet — ask us if you need a copy.

How to use it

Make up a fresh solution and renew it every 24 hours.

• Sterilising equipment & utensils

- 1 Wash items thoroughly in warm, soapy water first, then rinse with cold water.
- 2 Fill a non-metallic container with 5L of cold water and add **30ml (one capful)** of Milton fluid.
- 3 Fully submerge the items, making sure there are no trapped air bubbles.
- 4 Leave for 15 minutes. No need to rinse — items are ready to use straight away.
- 5 Items can stay in the solution until needed. Renew the solution every 24 hours.

• Disinfecting surfaces

- 1 Fill a container with 5L of cold water and add **90ml (three capfuls)** of Milton fluid.
- 2 Make sure the surface is clean first, then apply with a cloth, fully wetting the surface.
- 3 Leave for 5–15 minutes, then rinse — no need to rinse surfaces in direct food contact.
- 4 The solution stays effective for 24 hours.

Product details

Active ingredient	Sodium hypochlorite
Solution life	24 hours once made up

Suitable for: feeding equipment, utensils, non-metallic surfaces and items.
Not for metal surfaces.